



Rotavapor® R-220 EX / R-250 EX

EX-protection fulfilling highest safety levels

Intuitive operation and handling guarantee safe and efficient industrial distillation. Customize your Industrial Rotavapor® to your needs. Meeting the newest EX regulations, you are always on the safe side for distillation with a flask size of up to 50 liters.

Safe

According to ATEX regulation
for maximum protection



Easy

Intuitive and easy operation
at the highest safety level



Long lasting

Proven technology and
quality enhances lifetime

R-220 EX / R-250 EX: Your most important benefits



Safe

- ATEX approved for highest safety level
- Immediate shut-down due to safety cut-off
- Extend functionality and safety with a wide range of accessories
- Optional Ex-proof vacuum pump
- Designed for safe operation in harsh environment
- Safe handling and spring-loaded valves protect the glass



Easy

- Intuitive operation of the Rotavapor® and its integrated vacuum controller
- Excellent handling by the easy readable display
- Easy flask handling due to the unique flask coupling



Long lasting

- Long life time due to exclusive use of corrosion-free stainless steel
- Low maintenance by using highly chemical resistant materials
- Trouble-free operation by robust industrial valves

Complete your portfolio



Vacuum Pump
V-710

The powerful vacuum
source



Rotavapor®
R-220 SE

Best choice for
your scale up lab



Rotavapor®
R-215

Excellence in
evaporation



Recirculating Chiller
F-125

“Support center” for the
Rotavapor® R-220 SE

